

Monomethyl-Histone H3-K27 (H3K27Me1) Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Ihc, Chip & Elisa

Item Number: CM0000406



Introduction

Rabbit polyclonal antibody against monomethylated histone H3 at lysine 27 (H3K27me1). Validated for western blot, immunofluorescence, immunohistochemistry, dot blot, ChIP, and ELISA. Predicted cross-reactivity with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K27 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; MonoMethyl-Histone H3-K27
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic monomethylated peptide around K27 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, IHC-P, DB, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag