

Monomethyl-Histone H3-K27 Rabbit Mab (Q16695/P68431) - Methylated Histone H3 Antibody For Wb, If/Icc, Elisa

Item Number: CM0000255



Introduction

Rabbit monoclonal antibody specific to monomethylated histone H3 lysine 27 (H3K27me1) validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat, and plant histone H3, ideal for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K27 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; MonoMethyl-Histone H3-K27
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic monomethylated peptide around K27 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Arabidopsis thaliana, Rice, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag