


Monomethyl-Histone H3-K36 (Q16695 / P68431) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa, Chip-Seq

Item Number: CM0000316



Introduction

MonoMethyl-Histone H3-K36 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA, and ChIP-seq applications. Reacts with human, mouse, rat, and other wide range predicted species. Detects monomethylation at lysine 36 of histone H3, a key epigenetic mark in chromatin regulation.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K36 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; MonoMethyl-Histone H3-K36
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	A synthetic monomethylated peptide around K36 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag