

Monomethyl-Histone H3-K36 Rabbit Mab (Q16695/P68431) For Wb, Db, Ihc-P, Elisa

Item Number: CM0000386



Introduction

Rabbit monoclonal antibody detecting monomethylated lysine 36 of human histone H3 (SWISS Q16695/P68431). Predicted to cross-react with mouse, rat, and other species. Suitable for Western blot, dot blot, IHC-P, and ELISA.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K36 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; MonoMethyl-Histone H3-K36
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	A synthetic monomethylated peptide around K36 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag