

Monomethyl-Histone H3-K36 Rabbit Mab (Q16695/P68431) For Wb, Db, Ip, Elisa, Chip

Item Number: CM0000361



Introduction

MonoMethyl-Histone H3-K36 Rabbit mAb (SWISS: Q16695/P68431) for epigenetics research. Validated in WB, DB, IP, ELISA, ChIP; human, mouse, rat reactivity. 15kDa, ice-bag shipping.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K36 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; MonoMethyl-Histone H3-K36
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide A synthetic monomethylated peptide around K36 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag