

Monomethyl-Histone H3-K36 Rabbit Polyclonal Antibody (Q16695/P68431) - Wb, Chip, Elisa

Item Number: CM0000261



Introduction

Rabbit polyclonal antibody detecting mono-methylated histone H3 at lysine 36 (K36) for Western blot, IF/ICC, ELISA, ChIP, and ChIP-seq. Cross-reacts with human, mouse, rat, and predicted wide range species. Ideal for epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K36 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; MonoMethyl-Histone H3-K36
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag