


Monomethyl-Histone H3-K4 Rabbit Pab (Q16695/P68431) — Validated For Wb, Ihc-P, Chip & More

Item Number: CM0000364



Introduction

Rabbit polyclonal antibody specific to monomethylated histone H3 at lysine 4 (H3K4me1). Validated for WB, IHC-P, IF/ICC, ELISA, ChIP, and ChIP-seq. Predicted to cross-react across human, mouse, rat, and other species.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K4 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag