

Monomethyl-Histone H3-K79 Rabbit Polyclonal Antibody (Q16695 / P68431) For Wb, Elisa, Chip

Item Number: CM0000268



Introduction

Rabbit polyclonal antibody specific for monomethylated lysine 79 of histone H3 (SWISS Q16695/P68431). Validated for WB, ELISA, ChIP; cross-reacts with human, mouse, rat, and predicted wide range.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K79 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; MonoMethyl-Histone H3-K79
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide A synthetic monomethylated peptide around K79 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag