

Monomethyl-Histone H3-K79 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Ihc, Chip

Item Number: CM0000345



Introduction

Rabbit polyclonal antibody targeting monomethylated histone H3 at lysine 79 (K79me1), validated for WB, IHC-P, IF/ICC, ELISA, ChIP, and ChIP-seq. Cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K79 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; MonoMethyl-Histone H3-K79
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	A synthetic monomethylated peptide around K79 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag