


Monomethyl-Histone H3-K9 (Q16695/P68431) Rabbit Pab For Wb, Chip-Seq, Elisa

Item Number: CM0000375



Introduction

Rabbit polyclonal antibody specific for monomethylated histone H3 at lysine 9 (H3K9me1). Suitable for WB, ChIP-seq, IHC, ELISA. Cross-reacts with human, mouse, rat. Ideal for epigenetics and IVD raw material supply.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K9 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; MonoMethyl-Histone H3-K9
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide. A synthetic monomethylated peptide around K9 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag