


Monomethyl-Histone H3-K9 Rabbit Mab (Q16695/P68431) — Monoclonal Antibody For Epigenetic Research

Item Number: CM0000319



Introduction

MonoMethyl-Histone H3-K9 Rabbit mAb (Q16695/P68431) is a high-specificity monoclonal antibody for H3K9me1 detection in WB, DB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat, and other species predicted. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-K9 Rabbit mAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; MonoMethyl-Histone H3-K9
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide A synthetic monomethylated peptide around K9 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag