


Monomethyl-Histone H3-R17 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Ihc-P, If/Icc, Elisa

Item Number: CM0000409



Introduction

Rabbit polyclonal antibody targeting monomethylated R17 of human histone H3 (Q16695/P68431). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with Human, Mouse, Rat, and other species.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-R17 Rabbit pAb
Remarks/Alias	H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3Fj; H3C10; H3C11; HIST1H3j; MonoMethyl-Histone H3-R17
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag