


Monomethyl-Histone H3-R2 Rabbit Polyclonal Antibody (Q16695/P68431) | Wb, If/Icc, Elisa

Item Number: CM0000265



Introduction

Rabbit polyclonal antibody specific for monomethylated histone H3 at arginine 2. Suitable for WB, IF/ICC, ELISA with broad species cross-reactivity (human, mouse, rat). Ideal for epigenetics research and chromatin dynamics.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-R2 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; MonoMethyl-Histone H3-R2
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	A synthetic monomethylated peptide around R2 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag