


Monomethyl-Histone H3-R26 Rabbit Pab - Q16695/P68431 For Wb, If/Icc, Elisa

Item Number: CM0000389



Introduction

Rabbit polyclonal antibody specifically recognizing monomethylated Histone H3 at arginine 26 (Q16695/P68431). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat, and predicted wide-range species.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-R26 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; MonoMethyl-Histone H3-R26
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide / A synthetic monomethylated peptide around R26 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag