

Monomethyl-Histone H3-R8 (Q16695/P68431) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa

Item Number: CM0000404



Introduction

High-quality MonoMethyl-Histone H3-R8 rabbit polyclonal antibody suitable for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat histone H3.1t and H3.1 monomethylated at arginine 8, ~15 kDa.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H3-R8 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; MonoMethyl-Histone H3-R8
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	A synthetic monomethylated peptide around R8 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag