

Monomethyl-Histone H4-R3 Rabbit Pab (P62805) | If/Icc, Elisa Ready Antibody

Item Number: CM0000311



Introduction

Rabbit polyclonal antibody against monomethylated arginine 3 of human Histone H4 (SWISS P62805). Suitable for IF/ICC and ELISA; cross-reacts with human, mouse, rat, and other species. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H4-R3 Rabbit pAb
Remarks/Alias	H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4FA; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; H4C16; HIST1H4A; MonoMethyl-Histone H4-R3
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	A synthetic monomethylated peptide around R3 of human Histone H4 (NP_003529.1)
Source	Rabbit
Category	Methylated Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag