


Monomethyl-Histone H4-R3 Rabbit Polyclonal Antibody | P62805 | Wb, Ihc-P, If/Icc, Elisa

Item Number: CM0000403



Introduction

High-specificity rabbit polyclonal antibody detecting monomethylated arginine 3 of human histone H4 (SWISS P62805). Validated for Western blot, immunohistochemistry, IF/ICC, and ELISA; cross-reacts with human, mouse, rat, and predicted wide-range species. Supporting epigenetic research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-Histone H4-R3 Rabbit pAb
Remarks/Alias	H4/I; H4C1; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4FI; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; H4C16; HIST1H4B; MonoMethyl-Histone H4-R3
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	Synthetic peptide A synthetic monomethylated peptide around R3 of human Histone H4 (NP_003529.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag