

Monomethyl-Nono/P54Nrb (R357/R365) Rabbit Polyclonal Antibody (Q15233) For Wb, If/Icc, Elisa

Item Number: CM0000368



Introduction

Rabbit polyclonal antibody against monomethylated R357/R365 of human NONO/p54nrb (UniProt Q15233). Validated for WB, IF/ICC, ELISA, with cross-reactivity to human and rat. Ideal for studying NONO's roles in splicing, paraspeckles, DNA repair, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	MonoMethyl-NONO/p54nrb-R357/R365 Rabbit pAb
Remarks/Alias	P54; NMT55; NRB54; MRXS34; P54NRB; PPP1R114; MonoMethyl-NONO/p54nrb-R357/R365
Species	Human
GeneID (human)	4841
GeneID	4841
Immunogen	Synthetic peptide A synthetic monomethylated peptide around R357 & R365 of human NONO/p54nrb (Q15233).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q15233
Protein Weight	54kDa
Shipping	Ice bag