


Monomethyl-Yap1-K342 Rabbit Polyclonal Antibody (P46937) For Wb And Elisa

Item Number: CM0000259



Introduction

High-quality rabbit polyclonal antibody specifically recognizing monomethylated YAP1 at lysine 342. Ideal for WB and ELISA applications in human, mouse, and rat samples. Insights into Hippo signaling.

[Learn More](#)

Parameter	Description
Product Name	MonoMethyl-YAP1-K342 Rabbit pAb
Remarks/Alias	YAP; YKI; COB1; YAP2; YAP-1; YAP65; MonoMethyl-YAP1-K342
Species	Human
GeneID (human)	10413
GeneID	10413
Immunogen	A synthetic monomethylated peptide around K342 of human YAP1 (NP_001123617.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46937
Protein Weight	54kDa
Shipping	Ice bag