


Morf4L1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ubu8

Item Number: CM0026940



Introduction

High-quality rabbit polyclonal antibody targeting MORF4L1 (MRG15). Validated for WB and ELISA with human, mouse, rat reactivity. Suitable for chromatin regulation and DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	MORF4L1 Rabbit pAb
Remarks/Alias	Eaf3; MEAF3; MRG15; FWP006; S863-6; HsT17725; MORFRG15; MORF4L1
Species	Human
GeneID	10933
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human MORF4L1 (NP_006782.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBU8
Protein Weight	41kDa
Shipping	Ice bag