

Mouse Anti-Human Cd11B Monoclonal Antibody For Flow Cytometry - P11215

Item Number: CM0012060



Introduction

-conjugated mouse anti-human CD11b monoclonal antibody validated for flow cytometry. Targets integrin alpha-M (ITGAM), a leukocyte adhesion receptor in innate immunity. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD11b mAb
Remarks/Alias	CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6
Species	Human
GeneID	3684
Immunogen	A synthesized peptide derived from human CD11b
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P11215
Protein Weight	127kDa
Shipping	Ice bag