

Mouse Anti-Human Cd11B Monoclonal Antibody For Flow Cytometry - P11215

Item Number: CM0011932



Introduction

-conjugated mouse monoclonal antibody targeting human CD11b (integrin alpha-M). Suitable for flow cytometry applications, cross-reacts with human. Ideal for immunology research on monocyte/macrophage adhesion and phagocytosis.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD11b mAb
Remarks/Alias	CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6
Species	Human
Gene ID	3684
Immunogen	Synthetic peptide derived from human CD11b
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P11215
Protein Weight	127kDa
Shipping	Ice bag