

Mouse Anti-Human Cd11C Monoclonal Antibody For Fc - P20702

Item Number: CM0011473



Introduction

Mouse monoclonal antibody targeting human CD11c (integrin alpha-X), conjugated to . Validated for flow cytometry (FC), ideal for studying monocyte and granulocyte adhesion and inflammatory responses.

[Learn More](#)

Parameter	Specification
Product Name	Mouse anti-Human CD11c mAb
Remarks/Alias	CD11C; SLEB6
Species	Human
GeneID	3687
Immunogen	Recombinant protein of human CD11c.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20702
Protein Weight	128kDa
Shipping	Ice bag