

Mouse Anti-Human Cd15/Ssea-1 Monoclonal Antibody For Flow Cytometry - P22083

Item Number: CM0009843

Introduction

-conjugated mouse monoclonal antibody targeting human CD15/SSEA-1 for flow cytometry. Recognizes the Lewis X glycan epitope on glycoproteins; useful for studying stem cell and myeloid cell surface marker expression.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD15/SSEA-1 mAb
Remarks/Alias	LeX; CD15; ELFT; FCT3A; FUTIV; SSEA-1; FUC-TIV
Species	Human
GeneID (Human)	2526
GeneID	2526
Immunogen	A synthesized peptide derived from human CD15/SSEA-1.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P22083
Protein Weight	59 kDa
Shipping	Ice bag