

Mouse Anti-Human Cd15/Ssea-1 Monoclonal Antibody For Flow Cytometry - P22083

Item Number: CM0009496



Introduction

-conjugated mouse monoclonal antibody targeting human CD15/SSEA-1 (FUT4) for flow cytometric analysis. Validated for human samples, recognizes Lewis X antigen, key in leukocyte biology and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD15/SSEA-1 mAb
Remarks/Alias	LeX; CD15; ELFT; FCT3A; FUTIV; SSEA-1; FUC-TIV
Species	Human
Gene ID	2526
Immunogen	A synthesized peptide derived from human CD15/SSEA-1
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P22083
Protein Weight	59kDa
Shipping	Ice bag