

Mouse Anti-Human Cd3 Monoclonal Antibody For Fc And In Vitro T Cell Stimulation - P07766

Item Number: CM0009312



Introduction

Mouse monoclonal antibody targeting human CD3 epsilon chain. Validated for flow cytometry and in vitro T cell stimulation/activation. Reacts with human CD3. Suitable for immunology and T cell research.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon; CD3E Antigen
Species	Human
GeneID (Human)	916
GeneID	916
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC, in vitro T cell stimulation/activation
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag