

Mouse Anti-Human Cd3 Monoclonal Antibody For Flow Cytometry - P07766

Item Number: CM0010517



Introduction

Mouse anti-Human CD3 monoclonal antibody for flow cytometry, targeting the CD3 epsilon chain (P07766). Provides reliable detection of T lymphocytes in diagnostic and research settings.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
Gene ID (Human)	916
Immunogen	Recombinant protein of human CD3
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag