

Mouse Anti-Human Cd56 (Ncam1) Monoclonal Antibody For Flow Cytometry - P13591

Item Number: CM0011020



Introduction

Mouse monoclonal antibody against human CD56 (NCAM1), validated for flow cytometry. Ideal for immunophenotyping, stem cell research, and neuroscience applications. Detects cell adhesion molecule involved in neural development and virus entry.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD56 mAb
Remarks/Alias	CD56; NCAM; MSK39
Species	Human
GeneID (Human)	4684
GeneID	4684
Immunogen	Recombinant protein of human CD56
Source	Mouse
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human
SWISS	P13591
Protein Weight	40kDa/73kDa/80kDa/83kDa/93kDa/95kDa
Shipping	Ice bag