


Mouse Anti-Rabbit IgG (Fc) Monoclonal Antibody For Fc - P01870

Item Number: CM0004912



Introduction

High-quality -conjugated mouse monoclonal antibody targeting the Fc region of rabbit IgG (SWISS P01870), validated for flow cytometry. Ideal for immunoassays and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Rabbit IgG (Fc) mAb
Species	Rabbit
GeneID (human)	3500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 96-323 of rabbit IgG (Fc) (P01870)
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rabbit
SWISS	P01870
Protein Weight	35kDa
Shipping	Ice bag