


Mouse Anti-Rabbit IgG (Fc) Monoclonal Antibody For Flow Cytometry - P01870

Item Number: CM0010034



Introduction

-conjugated mouse monoclonal antibody targeting rabbit IgG Fc region (UniProt P01870). Ideal for flow cytometry, ensuring specific detection with recombinant protein-validated immunogen. High quality for immunoassay development.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Rabbit IgG (Fc) mAb
Species	Rabbit
GeneID (Human)	3500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 96-323 of rabbit IgG (Fc) (P01870)
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rabbit
SWISS	P01870
Protein Weight	35kDa
Shipping	Ice bag