


Mouse Anti-Rabbit IgG (Fc) Monoclonal Antibody For Flow Cytometry - P01870

Item Number: CM0006837



Introduction

High-specificity -conjugated mouse monoclonal antibody targeting rabbit IgG Fc region (P01870). Validated for flow cytometry; ideal for detecting rabbit immunoglobulins in immunoassays and immunofluorescence-based applications.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Rabbit IgG (Fc) mAb
Species	Rabbit
Gene ID (Human)	3500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 96-323 of rabbit IgG (Fc) (P01870)
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rabbit
SWISS	P01870
Protein Weight	35kDa
Shipping	Ice bag