

Mouse Anti-Rabbit Igm (Fc) Monoclonal Antibody For Fc - P03988

Item Number: CM0010952

Introduction

Mouse monoclonal antibody against rabbit IgM Fc region, conjugated to for flow cytometry. Targets the constant region of rabbit immunoglobulin mu heavy chain (P03988). Reacts specifically with rabbit.

[Learn More](#)

Parameter	Value
Product Name	Mouse Anti-Rabbit IgM (Fc) mAb
Species	Rabbit
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 107-458 of rabbit IgM (Fc) (P03988)
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rabbit
SWISS	P03988
Protein Weight	52kDa
Shipping	Ice bag