

Mouse Anti-Rabbit Igm (Fc) Monoclonal Antibody For Flow Cytometry - P03988

Item Number: CM0004833

Introduction

Mouse monoclonal -conjugated antibody specific for rabbit IgM Fc region (P03988). Suitable for flow cytometry, no cross-reactivity with other species. Ideal for rabbit immunoglobulin detection in research and diagnostic assays.

[Learn More](#)

Attribute	Value
Product Name	Mouse Anti-Rabbit IgM (Fc) mAb
Species	Rabbit
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 107-458 of rabbit IgM (Fc) (P03988)
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rabbit
SWISS	P03988
Protein Weight	52kDa
Shipping	Ice bag