

Mouse IgG 1 Isotype Control - Isotype Control

Item Number: CM0009303

Introduction

The Mouse IgG 1 isotype control is a monoclonal antibody conjugated to for use as a negative control in flow cytometry experiments. Sourced from mouse, it ensures accurate gating and background assessment.

[Learn More](#)

Parameter	Value
Product Name	Mouse IgG 1 isotype control
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag