

Mouse IgG 1 Isotype Control

Item Number: CM0005413

Introduction

Mouse IgG 1 isotype control, a monoclonal antibody for flow cytometry (FC). Ideal for establishing background staining in FC assays, sourced from mouse with recombinant immunogen. Shipped on ice bag.

[Learn More](#)

Parameter	Value
Product Name	Mouse IgG 1 isotype control
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag