

Mouse Igg 2A Isotype Control -

Item Number: CM0011618

Introduction

Mouse IgG 2a isotype control for flow cytometry (FC) applications. Recombinant protein immunogen, mouse monoclonal antibody. Ideal for assessing background staining and nonspecific binding in immunofluorescence studies.

[Learn More](#)

Parameter	Value
Product Name	Mouse IgG 2a isotype control
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag