

Mouse Igg 2A Isotype Control

Item Number: CM0012522

Introduction

Mouse IgG 2a isotype control is a monoclonal antibody derived from mouse, designed for use as an isotype control in flow cytometry (FC). It helps distinguish specific from non-specific binding, ensuring accurate experimental results. Shipped on ice bags.

[Learn More](#)

Parameter	Specification
Product Name	Mouse IgG 2a isotype control
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag