

Mouse Igg 2A Isotype Control - Isotype Control

Item Number: CM0010304

Introduction

The Mouse IgG 2a isotype control is a monoclonal antibody from mouse, designed as a negative control for flow cytometry (FC). It helps ensure accurate gating and eliminate background staining in research assays.

[Learn More](#)

Parameter	Value
Product Name	Mouse IgG 2a isotype control
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag