

Mouse Igg 2A Isotype Control

Item Number: CM0009392

Introduction

The Mouse IgG 2a isotype control is a monoclonal antibody conjugate designed for flow cytometry (FC) applications. It serves as a negative control to assess non-specific binding and background staining in mouse IgG 2a-based assays.

[Learn More](#)

Parameter	Value
Product Name	Mouse IgG 2a isotype control
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag