

Mouse IgG 2B Isotype Control Antibody - A25491

Item Number: CM0008090

Introduction

Mouse monoclonal IgG2b isotype control for flow cytometry. Validated for mouse samples, this antibody assesses non-specific binding in immunoassays and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	Mouse IgG 2b isotype control
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag