

Mouse IgG 2C Isotype Control - A26280

Item Number: CM0005551

Introduction

-conjugated mouse monoclonal IgG2c isotype control for flow cytometry (FC). Used to assess non-specific binding and background. Recombinant protein immunogen, mouse source. Ships on ice. Validate in your system.

[Learn More](#)

Parameter	Value
Product Name	Mouse IgG 2c isotype control
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag