

Mouse Igg1 Isotype Control

Item Number: CM0011627

Introduction

-conjugated Mouse IgG1 isotype control monoclonal antibody for flow cytometry (FC) applications. Produced in mouse, immunogen recombinant protein. Ideal for establishing background staining levels in IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	Mouse IgG 1 isotype control
Immunogen	Recombinant protein
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag