

Mouse Mdsc Panel Monoclonal Antibody - Mouse Mdsc Panel

Item Number: CM0006243

Introduction

Mouse MDSC Panel is a monoclonal antibody cocktail for flow cytometry analysis of mouse myeloid-derived suppressor cells. It provides specific reactivity for mouse samples, supporting immunology research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Mouse MDSC Panel
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag