

Mouse T Cells Panel Monoclonal Antibody For Flow Cytometry - Mouse T Cells Panel

Item Number: CM0004460

Introduction

The Mouse T Cells Panel monoclonal antibody is designed for flow cytometry (FC) applications, specifically targeting mouse T cells. This product is ideal for immunophenotyping and research on mouse T cell subsets.

[Learn More](#)

Parameter	Specification
Product Name	Mouse T Cells Panel
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag