

Mouse Treg Panel Monoclonal Antibodies For Flow Cytometry - Mouse Treg Panel

Item Number: CM0012180

Introduction

High-quality mouse Treg panel monoclonal antibodies from rabbit for flow cytometry. Ideal for mouse T regulatory cell detection and analysis. Reliable IVD raw material from CamelBio.

[Learn More](#)

Parameter	Specification
Product Name	Mouse Treg Panel
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag