

Mouse Treg Panel Monoclonal Antibody (Rabbit) For Flow Cytometry - Mouse Treg Panel

Item Number: CM0008509

Introduction

Mouse Treg panel monoclonal antibody (rabbit) specifically for flow cytometry (FC) analysis of mouse regulatory T cells. Ideal for immunology research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Mouse Treg Panel
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag