

Mouse Ubn2 Rabbit Polyclonal Antibody - Ubn2

Item Number: CM0016970

Introduction

Rabbit polyclonal antibody targeting mouse UBN2 (NP_796159.3), raised against a synthetic peptide within amino acids 1-100. Suitable for Western blot and ELISA, with reported cross-reactivity to human UBN2. Predicted molecular weights range from 130 to 141 kDa.

[Learn More](#)

Parameter	Value
Product Name	UBN2 Rabbit pAb
Remarks/Alias	mKIAA2030, 2900060J04Rik, 6030408G03Rik, D130059P03Rik, UBN2
Species	Mouse
Gene ID (Human)	254048
Gene ID (Mouse)	320538
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse UBN2 (NP_796159.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
Protein Weight	130kDa/131kDa/141kDa
Shipping	Ice bag