

Mouse Wdr33 Rabbit Polyclonal Antibody - Wdr33

Item Number: CM0012912

Introduction

Rabbit polyclonal antibody targeting mouse WDR33, validated for WB and ELISA. Cross-reacts with human and rat WDR33, detecting isoforms at 30, 38, and 145 kDa. Ideal for protein research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	WDR33 Rabbit pAb
Remarks/Alias	WDC146; 1110001N06Rik; 2310011G05Rik; 2810021O11Rik; 8430413N20Rik; WDR33
Species	Mouse
Gene ID (Human)	55339
Gene ID	74320
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of mouse WDR33 (NP_001164437.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
Protein Weight	30kDa/38kDa/145kDa
Shipping	Ice bag