


Mpc1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Y5U8

Item Number: CM0023902



Introduction

MPC1 rabbit polyclonal antibody validated for WB, IP, ELISA with human, mouse, rat cross-reactivity. Targets mitochondrial pyruvate carrier 1, mediating pyruvate uptake and essential for cellular metabolism.

[Learn More](#)

Parameter	Value
Product Name	MPC1 Rabbit pAb
Remarks/Alias	MPYCD; BRP44L; CGI-129; SLC54A1; MPC1
Species	Human
Gene ID (Human)	51660
Gene ID	51660
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-109 of human MPC1 (NP_057182.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5U8
Protein Weight	12kDa
Shipping	Ice bag