

Mphosph6 Rabbit Polyclonal Antibody For Wb, Elisa - Q99547

Item Number: CM0014413



Introduction

Rabbit polyclonal antibody against human MPHOSPH6 (MPP6), validated for WB and ELISA. Recognizes human, mouse, and rat MPHOSPH6. Useful for studying RNA exosome function and rRNA processing.

[Learn More](#)

Parameter	Value
Product Name	MPHOSPH6 Rabbit pAb
Remarks/Alias	MPP; MPP6; MPP-6; MPHOSPH6
Species	Human
GeneID (Human)	10200
GeneID	10200
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human MPHOSPH6 (NP_005783.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99547
Protein Weight	19kDa
Shipping	Ice bag